

## Product datasheet

## Mouse Anti DXd Antibody[Biotin], mAb

Cat: NA0104

Name: Mouse Anti DXd Antibody[Biotin], mAb

Immunogen: DXd conjugated to keyhole limpet haemocyanin(KLH)

Conjugate: Biotin

Host Species: Mouse

Specificity: This antibody recognizes DXd and DXd structurally similar molecules. The antibody can bind to DS-8201.

Form: Liquid

Concentration: 0.5 mg/ml

Storage Buffer: 1×PBS, pH7.4, containing 1% BSA, 50% glycerol and 0.03% Proclin 300.

Clonality: Monoclonal

Isotype: Mouse IgG1,λ

Clone ID: 10G8

Shipping and Storage: Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.

Expiration Date: Specified on product label.

## Applications

Application	Recommended Usage
ELISA	0.1-2 µg/ml
Other applications	User-optimized

Note: Ningyang Biotech can offer this product according to your requirement, including product size, buffer components, etc.

**For research use only. Not intended for human and animal therapeutic or diagnostic use.**