

## Product datasheet

## Mouse Anti DYKDDDDK Tag Antibody, mAb

Cat: NA0113

Name: Mouse Anti DYKDDDDK Tag Antibody, mAb

Immunogen: DYKDDDDK peptide conjugated to keyhole limpet haemocyanin(KLH)

Host Species: Mouse

Specificity: This antibody recognizes C-terminal, N-terminal, and internal region of DYKDDDDK tag fusion proteins.

Form: Liquid

Concentration: 0.5 mg/ml

Storage Buffer: 1×PBS, pH7.4, containing 0.1% BSA, 50% glycerol and 0.03% Proclin 300.

Conjugate: Horseradish peroxidase (HRP)

Clonality: Monoclonal

Isotype: Mouse IgG2b,  $\kappa$

Clone ID: 66F10

Shipping and Storage: Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.

Expiration Date: Specified on product label.

## Applications

Application	Recommended Usage
ELISA	0.1-1 $\mu$ g/ml
Western Blot	0.5-1 $\mu$ g/ml
Other applications	User-optimized

Note: Ningyang Biotech can offer this product according to your requirement, including product size, buffer components, etc.

**For research use only. Not intended for human and animal therapeutic or diagnostic use.**