

Product datasheet

Mouse Anti His Tag Antibody[HRP], mAb

Cat: NA0106

Name: Mouse Anti His Tag Antibody[HRP], mAb

Immunogen: HHHHHH peptide conjugated to keyhole limpet haemocyanin(KLH)

Conjugate: Horseradish peroxidase (HRP)

Host Species: Mouse

Specificity: This antibody recognizes C-terminal, N-terminal, and internal region of His tag fusion proteins.

Form: Liquid

Concentration: 0.5 mg/ml

Storage Buffer: 1×PBS, pH7.4, containing 1% BSA, 50% glycerol and 0.03% Proclin 300.

Clonality: Monoclonal

Isotype: Mouse IgG1, κ

Clone ID: 90B8

Shipping and Storage: Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.

Expiration Date: Specified on product label.

Applications

Application	Recommended Usage
ELISA	0.05-0.5 μ g/ml (1:1000-1:10,000)
Western Blot	0.1-1 μ g/ml (1:500-1:5,000)
Other applications	User-optimized

Note: We can offer this product according to your requirement, including product size, buffer components, etc.

For research use only. Not intended for human and animal therapeutic or diagnostic use.