

Product datasheet

Rabbit Anti-Camelid VHH[HRP],pAb

Cat: NA0024

Name: Rabbit Anti-Camelid VHH[HRP],pAb

Immunogen: Camelid VHH

Conjugate: Horseradish peroxidase (HRP)

Host Species: Rabbit

Specificity: This antibody recognizes camelid VHH, recombinant VHH and humanized VHH.

Form: Liquid

Concentration: 0.5 mg/ml

Storage Buffer: 1×PBS, pH7.4, containing 1% BSA 50% glycerol and 0.03% Proclin 300.

Clonality: Polyclonal

Isotype: Rabbit IgG

Shipping and Storage: Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.

Expiration Date: Specified on product label.

Applications

Application	Recommended Usage
ELISA	1:5,000-1:50,000
Western Blot	1:1,000-1:10,000
Other applications	User-optimized

Note: Ningyang Biotech can offer this product according to your requirement, including product size, buffer components, etc.

For research use only. Not intended for human and animal therapeutic or diagnostic use.