

## Product datasheet

## Rabbit Anti E Tag Antibody, pAb

Cat: NA0067

Name: Rabbit Anti E Tag Antibody, pAb

Immunogen: E tag peptide(GAPVPYPDPLEPR) conjugated to keyhole limpet haemocyanin(KLH)

Host Species: Rabbit

Specificity: This antibody recognizes C-terminal and N-termina region of E Tag fusion proteins.

Form: Liquid

Concentration: 0.5 mg/ml

Storage Buffer: 1×PBS, pH7.4, containing 50% glycerol and 0.03% Proclin 300.

Purification: Antigen affinity purification

Clonality: Polyclonal

Isotype: Rabbit IgG


Shipping and Storage: Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.

Expiration Date: Specified on product label.

## Applications

Application	Recommended Usage
ELISA	0.2-1µg/ml
Western Blot	0.5-1µg/ml
Other applications	User-optimized

Note: Ningyang Biotech can offer this product according to your requirement, including product size, buffer components, etc.



**Figure 1.** Detection of E tag in recombinant fusion protein were analyzed by Western blot using Rabbit Anti E Tag Antibody, pAb (NingyangBiotech,Cat.NA0067, at 1  $\mu$ g/mL concentration).

The signal was developed with Goat Anti-Rabbit IgG(H&L)[HRP], pAb(NingyangBiotech Cat.NA0030).

**For research use only. Not intended for human and animal therapeutic or diagnostic use.**