

Product datasheet

Rabbit Anti GST Tag Antibody[HRP], pAb

Cat: NA0057

Name: Rabbit Anti GST Tag Antibody[HRP], pAb

Immunogen: Recombinant Glutathione-S-Transferase(GST)

Conjugate: Horseradish peroxidase (HRP)

Host Species: Rabbit

Specificity: This antibody recognizes GST tag fusion proteins.

Form: Liquid

Concentration: 0.5 mg/ml

Storage Buffer: 1×PBS, pH7.4, containing 1%, 50% glycerol and 0.03% Proclin 300.

Clonality: Polyclonal

Isotype: Rabbit IgG

Shipping and Storage: Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.

Expiration Date: Specified on product label.

Applications

Application	Recommended Usage
ELISA	0.05-0.5µg/ml
Western Blot	0.2-1µg/ml
Other applications	User-optimized

Note: Ningyang Biotech can offer this product according to your requirement, including product size, buffer components, etc.




Figure 1. Detection of GST tag in recombinant fusion protein were analyzed by Western blot using Rabbit Anti GST Tag Antibody[HRP], pAb (NingyangBiotech,Cat.NA0057, at 1 μ g/mL concentration).

For research use only. Not intended for human and animal therapeutic or diagnostic use.