

Streptavidin-HRP

Cat: NP004

Name: Streptavidin-HRP

Size: 1 mg

Conjugate: Horseradish peroxidase (HRP)

Form: Liquid

Concentration: 1 mg/ml

Storage Buffer: 1×PBS, pH7.4, containing 1% BSA, 50% glycerol and 0.03% Proclin 300.

Shipping and Storage: Shipped in ice bag. Store in the dark. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.

Expiration Date: Specified on product label.

Applications

Application	Recommended Usage
ELISA	1:5,000-1:20,000
Western Blot/Dot Blot	1:2,000-1:20,000
Other applications	User-optimized

Note: Ningyang Biotech can offer this product according to your requirement, including product size, buffer components, etc.

For research use only. Not intended for human and animal therapeutic or diagnostic use.